



E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887


FAX: +886-2-2537-1870

EMAIL: service@echromtech.com

ADD: 10F-12, NO.20, Min Chuan W. Road, Taipei 104,
Taiwan

WEB: www.echromtech.com

A07- HPLC Filter Holder



FS-10, FS-20



FS-01

Descriptions

Widely used in gravimeter, microanalysis, colloid dispersion and bacteria-free experiments, this HPLC filter holder can filtrate and degas the mobile phase of HPLC analysis. Besides playing a critical role in lifetime extension of apparatus and HPLC column, the filter holder can enhance measurement accuracy.

The filtration bottles have high pressure endurance, fabulous interchange ability, and well-proportional wall. They are made according to international standard sizes so that they can match most products of foreign brands. In addition, the newly designed solvent pickup adaptor can filtrate or collect solvent directly from solvent bottle through film of filter. This new design makes the process of filtrating become easier and it can dramatically increase effectiveness and reduce air pollution.

Features

- A. High-quality and hard-glass materials that can endure more than 200 temperature difference.
- B. Elegantly made with well-proportional wall.
- C. Good endurance of high pressure and fabulous interchange ability.
- D. Made according to international standard sizes that can match most foreign brand products.

Specifications

Packing list of FS-10 and FS-20

Glass Funnel: 47mm, 300ml

Plug

Fritted Glass Support Base: 47mm

Aluminum Clamp: 47mm

Erlenmeyer Flask: 1000ml (FS-10), 2000ml (FS-20)

Packing list of FS-01

Glass Funnel: 47mm, 300ml

Silicone Plug

Aluminum Clamp: 47mm

Fritted Glass Support Base: 47mm

Glass Cover

Beaker

Teflon Base