



## E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887

FAX: +886-2-2537-1870

EMAIL: [service@echromtech.com](mailto:service@echromtech.com)

ADD: 10F-12, NO.20, Min Chuan W. Road,  
Taipei 104, Taiwan

WEB: [www.echromtech.com](http://www.echromtech.com)

### B06- Shaker



**SH-50AFS**



**SH-50ARCS**



**SH-100AVS**

### Descriptions

- A. There are various types of shakers are available, including orbital shaking motion, circular shaking motion, and horizontal shaking motion.
- B. Time adjustment: operation of set time (3 ~ 60min) or continuous.
- C. Shakers are made of materials of high thermostability, anti-acid and anti-alkali materials, having good endurance.
- D. Shaker base has four stationary suckers to prevent any migration of work.

### Features

- A. Adjustable speed from 0 to 150rpm.
- B. Able to be placed in a temperature chamber with max temperature to 80°C
- C. Especially ideal for fermentation and cell culture.



## E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887

FAX: +886-2-2537-1870

EMAIL: [service@echromtech.com](mailto:service@echromtech.com)

ADD: 10F-12, NO.20, Min Chuan W. Road,  
Taipei 104, Taiwan

WEB: [www.echromtech.com](http://www.echromtech.com)

### Specifications

Model	SH-50AFS	SH-50ARCS	SH-100AVS
Speed Range	0 ~ 150rpm	0 ~ 350rpm	0 ~ 2500rpm
Vibrating Motion	Orbital	Horizontal	Circular
Vibrating Displacement	--	10mm	4.5mm
Max Load	2Kg	3Kg	1Kg
Time Set	3~60min or continuous	3~60min or continuous	3~60min or continuous
Weight (Kg)	2.4Kg	2.7Kg	8.2Kg

Size of plate: 21cm\*29cm

Plate assembly: conical flask, 50ml \*8 sets, or 100ml \*5sets, or 200ml \*2 sets

Standard stainless holder with our products: 50ml \*2, 100ml \*2 and 200ml \*2

### The equipped with the holders (B) and (D)

