



E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887

FAX: +886-2-2537-1870

EMAIL: service@echromtech.com

ADD: 10F-12, NO.20, Min Chuan W. Road, Taipei
104, Taiwan

WEB: www.echromtech.com

B12- FA/JA Electronic Balance



Features

- A. Automatic failure detection function.
- B. Automatic zero setting function with brand-new outward appearance.
- C. STN wide visual angle back-lighted liquid crystal display (LCD) has superior performance.
- D. RS-232C output interface.

Specifications

Model	FA-1004	FA-1104	FA-1604	FA-2004	FA-2104
Capacity	100g	110g	160g	200g	210g
Minimum Display	0.0001g	0.0001g	0.0001g	0.0001g	0.0001g
Pan Size	Φ90mm	Φ90mm	Φ90mm	Φ90mm	Φ90mm

Model	JA-1103	JA-2103	JA-3103	JA-5103	JA-10003
Capacity	110g	210g	310g	510g	1000g
Minimum Display	0.001g	0.001g	0.001g	0.001g	0.001g
Pan Size	Φ90mm	Φ90mm	Φ90mm	Φ110mm	Φ110mm

Packing list

- A. Electronic balance
- B. Stainless steel pan
- C. Power line
- D. Wind-cover
- E. Wind-top
- F. User manual



E-Chrom Tech Co., Ltd.

TEL: +886-2-2511-1887

FAX: +886-2-2537-1870

EMAIL: service@echromtech.com

ADD: 10F-12, NO.20, Min Chuan W. Road, Taipei
104, Taiwan

WEB: www.echromtech.com

B12- SL Electronic Balance



Features

- A. Brand-new design with STN wide visual angle back-lighted liquid crystal display (LCD)
- B. Able to transfer to many measurement units.
- C. RS-232C output port as an optional equipment.

Specifications

Model	SL-102	SL-202	SL-302	SL-502	SL-602	SL-1002	SL-2002
Capacity	100g	200g	300g	500g	600g	1000g	2000g
Readability	0.01g	0.01g	0.01g	0.01g	0.01g	0.01g	0.01g
Pan Size	Φ110mm	Φ110mm	Φ110mm	Φ128mm	Φ128mm	Φ165mm	Φ165mm

Packing list

- A. Electronic balance
- B. Stainless steel pan
- C. Power line
- D. Wind-cover
- E. Wind-top
- F. User manual